

NEXUS™ Aortic Arch Stent Graft System

Endovascular repair of the aortic arch made easy and reliable



The only
**CE approved
branched device**
for endovascular
repair of the
aortic arch


NEXUS™ is manufactured by Endospan Ltd.



CryoLife®
Life Restoring Technologies®

Designed for the aortic arch

NEXUS™ Aortic Arch Stent Graft System is the innovative solution that transforms complex aortic arch surgical repair into a minimally invasive endovascular procedure.



Designed for long lasting effectiveness

Every detail of the NEXUS™ Aortic Arch Stent Graft System has been designed to provide durable effective results. Disconnection of the modules is prevented by the proprietary locking system, where the locking latches engage the inner stent graft of the docking sleeve.



Designed for safety

Since neurological complications are among the most devastating complications of TEVAR in the arch, NEXUS™ has been designed with many features to minimize this risk.

Lower profile 20 F system for greater access confidence

Can be introduced in vessels diameter as small as 7 mm reducing concerns for access complications

Anatomically pre-shaped delivery systems


Pre-shaped Main Module delivery system

Pre-curved Ascending Module delivery system

The main body is advanced using a through-and-through technique between the main access and the right subclavian artery, avoiding manipulation in the right common carotid artery

Dual High Flush ports for reduced embolic risks

Dual High Flush ports maximise clearance of residual air from delivery system




Simple and controlled 2-stage procedure


The implantation consists of a first surgical stage, where the left common carotid (and left subclavian) artery is rerouted, followed later by the endovascular stage.

Deployment steps


Main Module deployment




Ascending Module insertion & positioning



Ascending Module deployment



Final result



Ordering information

NEXUS™ Aortic Arch Stent Graft System is available in a range of sizes that can adapt to a wide range of anatomies.


Main Module*

	Diameter (mm)	Length (mm)
Branch	14	20
	17	30
	20	40
Main Body	32	180
	36	180
	40	180
	44	180

Ascending Module*

	Diameter (mm)	Length (mm)
Oriented	36	40
	40	55
	43	70

* all possible combinations of diameter and length are available



NEXUS™ Anatomical indications

Ascending Aorta


- Diameter of 29-39 mm
- Landing zone length of at least 30 mm

Descending Aorta

- Diameter of 26-40 mm
- Landing zone length of at least 30 mm

Brachiocephalic trunk

- Diameter of 11.5-18.5 mm
- Landing zone length of at least 20 mm
- Take off angle between the brachiocephalic artery and the aortic arch perpendicular should be $\geq 125^\circ$



NEXUS™ – the only CE marked endovascular branch system for the aortic arch. Ready when you are!

EU Authorized Representative:
MedNet GmbH, Borkstrasse 10, 48163 Münster, Germany

CE 0344

JT-BR-0050200-EN V01 12/2019

JOTEC GmbH, a fully owned subsidiary of CryoLife, Inc.

Lotzenäcker 23 • 72379 Hechingen, Germany
P +49 (0)7471 922-0 • F +49 (0)7471 922-100
info.europe@cryolife.com • www.jotec.com

© 2019 JOTEC GmbH, Germany. All rights reserved.

CryoLife, Inc. • 1655 Roberts Blvd., NW • Kennesaw, GA 30144 • USA
P +1 770 419-3355 • +1 800 438-8285 • F +1 770 590-3739
CustomerService@cryolife.com • www.cryolife.com



Endospine Ltd.
Maskit St. 4, Herzlia Business Park, Herzlia,
ISRAEL 46733



CryoLife®
Life Restoring Technologies®